

SurgiSLUSH™

Works FOR You

Free your team from making and maintaining slush during complex procedures.

Enjoy a cleaner,
quieter,
and more focused sterile field.

Get More from Your Slush System



Automated 4-Liter



Reusable
High-Tech

Protected
Sterility

°C Change
SURGICAL

www.cchangesurgical.com (877) 989-3737



Reusable Container inside Protective Container System™ with Locking, Tamper-Evident Cap.



With PCS™ You KNOW Your Slush is Sterile During Use.

Specifications:

- Specified C Change Surgical Containers only
- Power: 120 Volt A/C, 60 Hz, 15 Amps
- Main AC Power Fuse: 250VAC 15 amp time delay fuse (Class J)
- DC Power Fuse: 250VAC 5 amp fast-acting fuse (5 X 20 mm)
- Uses (approximately) 2,000 BTUs Per Hour

Ordering Information:

orders@cchangesurgical.com

CLINICAL REVIEWS & COMMENTS:

The slush was perfect!!!! (AND THE ROOM WAS SO QUIET) without that noisy slush making machine that we have been using! Thank you so much for your help and all of the information you provided us." BSN RN CNOR Unit Based Educator CVOR

This is the best slush I've ever used. Chair, Dept of Cardiothoracic Surgery, Academic Medical Center

Should slush in a sterile self-contained package be standard of care? I believe it should. Great patient-safety centered innovation! Medical Director and Practicing Urologist

After 30 years somebody finally figured out surgical slush. Cardiac Surgeon, Regional Medical Center

Nice job on the design. You guys are amazing. Transplant Nurse Manager - Academic Medical Center

Our experience has been very favorable. Chief Medical Officer, OPO Executive VP, Leading Transplant Surgeon

We have one...a fantastic device. OPO CEO

I love it. Recovery Manager - OPO

High Output. Auto-Controlled.

A SMART SYSTEM THAT WORKS FOR YOU. LITERALLY.

While your clinical teams get things done, smart SurgiSLUSH™ freezers automatically produce and maintain smooth slush in locked, protected containers. Smart freezers efficiently supply multiple procedures. Silent containers replace machines in the sterile field and save valuable space. Enjoy a cleaner, quieter and more focused sterile field. Get More from your slush system.

FAR LOWER EQUIPMENT SPEND - One Freezer Will Supply Multiple Rooms

MORE USEFUL - Secure, Portable Slush Supplies Procurements, Transplants, Cardiac, Urology, Hybrid Rooms - All Slush Applications. Consolidate!

EFFICIENT - Users Love the Auto-Start Option: Load Containers at the End of the Day. Auto-Start Automatically Produces Fresh, Ready-to-Use Slush Each Morning.

VERIFIABLE STERILITY AT TIME OF USE - Tamper-Evident, Hard-Shell Container System Locks in and Protects Sterile Slush

CLEAN - No Prolonged Slush Exposure to Ambient Environment; No Air Exchanging Machines in the Sterile Field

QUIET - Capped, Space-Saving Slush Containers Make NO Noise

REDUCE MEDICAL WASTE - And Save Sterile Fluid One Protected Liter at a Time.

SCALABLE - Add Containers to Increase Supply - Not Expensive Machines

COST-EFFECTIVE - Containers Replace Machines and Slash Costs; Locked Containers are Usable for 24 Hours and Re-slushed Easily

ALWAYS SMOOTH SLUSH - Reusable, Hydrophobic Containers Use Applied Materials Science to Dependably Produce Smooth Slush

www.cchangesurgical.com