

Programmable Auto-Freezers with Protected Containers

SurgiSLUSH™

Works FOR You

Free your team
from making and maintaining slush
during complex procedures.

Enjoy a more focused sterile field.

Protect Every Patient
Simplify Your Sterile Field
Automate Expertise



°C Change
SURGICAL

www.cchangesurgical.com (877) 989-3737



AutoSoft™ Reusable Slush Container inside ChillVault™ Protective Container with Locking, Tamper-Evident Cap.



With ChillVault™ Protective Container, You KNOW Your Slush is Sterile During Use.

Specifications:

- Specified C Change Surgical Containers only
- Power: 120 Volt A/C, 60 Hz, 15 Amps
- Main AC Power Fuse: 250VAC 15 amp time delay fuse (Class J)
- DC Power Fuse: 250VAC 5 amp fast-acting fuse (5 X 20 mm)
- Uses (approximately) 2,000 BTUs Per Hour

Ordering Information:

orders@cchangesurgical.com

CLINICAL REVIEWS & COMMENTS:

"The SurgiSLUSH™ freezer is a major improvement over our old system! This system can provide slush for multiple rooms. The OR environment is much quieter, and less crowded. We make the slush in the morning before our cases start so they are ready for the day and the quality of the slush is far superior over the other technology. With this system, our surgeons have smoother, high quality slush to work with and our nurses don't have to attend to making sure the slush is maintained ..." Clinical Coordinator | CVOR

"My team and our surgeons ... big supporters of SurgiSLUSH™!" Director | Surgical Services

"The chopping / pounding / mixing and timing issues of producing "acceptable sterile slush" comprise some of the more crude things going on in an OR today. It's the 21st century ... there is a better way. THIS!" CVOR Clinical Consultant

"Should slush in a sterile self-contained package be standard of care? I believe it should. Great patient-safety centered innovation!" Medical Director | Practicing Urologist

"Those impervious drapes get holes in them – can't tell you how many times I have had to tell a family about a potential exposure This is why I want a better system ... where sterility is never a concern." Chief, Cardiothoracic Surgery | VA Medical Center

"The slush was perfect!!!! (AND THE ROOM WAS SO QUIET)" BSN RN CNOR | Educator | CVOR

"This is the best slush I've ever used." Chair | Dept of Cardiothoracic Surgery | Academic Medical Center

The Best Advances Are Automatic.

Slush-using clinicians overwhelmingly requested clinical and economic improvements. We listened.

Teams embrace SurgiSLUSH™ because it automatically ensures patient safety, improves efficiency, increases utility and delivers significant savings.

- Closed system for infection prevention. No prolonged open exposure.
- Tamper-evident, sealed containers confirm sterility BEFORE each use.
- Atraumatic, auto-soft slush every time via hydrophobic containers.
- One touch to produce and maintain optimal slush.
- Smart freezers work FOR you while you focus on priorities.
- Highly efficient during complex cases. No tending slush units.
- Higher output. Slush and store. Make more.
- One freezer can supply multiple procedures with portable containers.
- Scalable. Increase supply with containers, not costly machines.
- Delay-start. Auto-produce slush at the time you choose for max efficiency.
- Sealed containers are usable for 24 hours.
- Slush efficiently pre-made outside of the sterile field.
- More sterile space. Sterile containers replace draped machines.
- Erase a trip hazard. One less cord in sterile field.
- Containers make no noise. More sterile quiet.
- Go green. Slash material waste with small containers.
- Lower equipment cost to increase accessibility.
- Consolidate purchasing. Supply ALL slush uses.
- Save sterile fluid using one liter at a time.

Clinician-Guided, Patient-Focused Design

Advance your infection prevention efforts automatically. Our closed, automated system ends troublesome sterile barrier failures and erases prolonged, ambient exposure. Confirm sterility before use, for every patient, with a modern system designed to meet today's standards.

Automate expert slush and simplify your sterile field. Outside of the sterile field, one-touch freezers automatically produce & maintain atraumatic slush in sealed, 24-hour containers. Efficient slush is ready when needed, allows you to focus on priorities and facilitates quiet focus. No more fussing with slush during complex cases. You also gain sterile space and eliminate a cord on the floor. Automated advances deliver easy wins.

Supply more slush with less equipment. Just slush, store and make more. Auto-produce slush at the time you choose for max efficiency. A closed system delivers important clinical, economic and practical advantages.

One-Touch. Protected. Auto-Soft. Hands-Free.

SurgiSLUSH™

Exceeding Today's Standards.
Automatically.

C Change
SURGICAL

www.cchangesurgical.com

Email info@cchangesurgical.com

Toll-free (877) 989-3737

HealthTrust | GSA | MAGNET Group



ChillVault™ Protective Container
AutoSoft™ Reusable Container

2-Liter and 4-Liter
SurgiSLUSH™ Programmable Auto-Freezer

Best-Practice Sterile Security Locks In Sterility. Automatically.



ChillVault™ Protective Containers

Locks & Protects Because You Need to Know

**ChillVault™ tamper-evident cap seals AutoSoft™ reusable slush container inside.
Confirm sterility BEFORE use ... for every patient, for every liter.
Automatically.**



www.cchangesurgical.com Email info@cchangesurgical.com Toll-free (877) 989-3737
HealthTrust | GSA | MAGNET Group